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

A screenshot of the 'MS/MS Ions Search' dialog box. The interface is filled with various input fields and dropdown menus, each marked with a red circle and a number from 1 to 14. 1: URL field (http://www.matrixscience.com/cgi/nph-mascot.exe?1). 2: User Name field (user). 3: Email field (user@bruker-daltonics.jp). 4: Search title field. 5: Taxonomy dropdown (All entries). 6: Database dropdown (NCBIInr). 7: Enzyme dropdown (Trypsin). 8: Fixed Modifications list (Acetyl (N-term), Amide (C-term), Biotin (K), Biotin (N-term), Carbamidomethyl (C)). 9: Variable Modifications list (AB\_old\_ICATd0 (C), AB\_old\_ICATd8 (C), Acetyl (K), Acetyl (N-term), Amide (C-term)). 10: Protein mass input field. 11: Missing cleavages max dropdown (1). 12: Peptide tol. ± input field (150 ppm). 13: MS/MS tol. ± input field (5 Da). 14: Instrument dropdown (Default). At the bottom, there are buttons for 'Copy Peaklist', 'Copy Masslist', 'Save as default', 'Start', and 'Exit'. A 'Peaklist' section contains a list of m/z values and their relative intensities.

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
**Mascot Search Results: Protein View - Microsoft Internet Explorer**

Hit to: [gi|130794280](#) Score: 310 Expect: 3.1e-27  
albumin [Bos taurus]

Nominal mass ( $M_n$ ): 71274; Calculated pI value: 5.82  
NCBI BLAST search of [gi|130794280](#) against nr  
Unformatted [sequence string](#) for pasting into other applications

Taxonomy: [Bos taurus](#)  
Links to retrieve other entries containing this sequence from NCBI Entrez:  
[gi|12190337](#) from [Bos taurus](#)  
[gi|13336842](#) from [Bos taurus](#)  
[gi|123307791](#) from [Bos taurus](#)

Fixed modifications: Carbamidomethyl (C)  
Cleavage by Trypsin: cuts C-term side of RR unless next residue is P  
Number of mass values searched: 56  
Number of mass values matched: 30  
Sequence Coverage: 53%

Matched peptides shown in **Bold Red**

```

1 MKVVTITSL LLFSSAYSRG VFRDTHSE IAHFKDLGE EHWKGLVIA
51 FSQYLQCCPF DEHKLVNEL IEFARTQVAD ESHAGCEKSL HILFGDELCK
101 VASLRETYGD NADCCERQEP ERNECFLSHK DSDPDLPKLK PDPNLTCDEF
151 KADEKFFVGR YLVEIARHP YFAPPELLY AHKYGQVFE CQAEKDGAC
201 LFKIETIHR FVLTSARQR LRCASIQRF ERALKAWVA RLSQRFPAE
251 FVVTIKLTDV LTKVHKCH GDVIFCDDK ADLAKYTCN GDTTSKIKR
    
```

Chordata: Craniata; Vertebrata; Euteleostomi; Cetartiodactyla; Ruminantia; Pecora; Bovidae;

Bovinae; Bos.

REFERENCE 1 (residues 1 to 607)  
AUTHORS Higer,C., Geigioni,F., De Beaufort,C., Michei,G., Freilinger,J. and Henzges,F.  
TITLE Differential binding of IgG and IgA antibodies to antigenic

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

A screenshot of the 'MS/MS Ions Search' dialog box. The window title is 'MS/MS Ions Search'. It contains several input fields and dropdown menus. The 'URL' field is set to 'http://Djkt951s/mascot/cgi/nph-mascot.exe?1'. The 'Database' dropdown is set to 'NCBI nr'. The 'Enzyme' dropdown is set to 'None'. The 'Fixed Modifications' list includes 'Acetyl (K)', 'Acetyl (N-term)', 'Amide (C-term)', 'Biotinylated (K)', and 'Biotinylated (N-term)'. The 'Variable Modifications' list includes 'Acetyl (K)', 'Acetyl (N-term)', 'Amide (C-term)', 'Biotinylated (K)', and 'Biotinylated (N-term)'. The 'Protein mass' field is empty, and the 'Peptide tol. ±' is set to '150 ppm'. The 'Charge state' is set to '1+', and the 'm/z' is set to '2061.830100'. The 'Data file' field is empty. The 'Peaklist' field contains a list of m/z values: 60.010613 205.720000, 86.251503 145.580000, 101.154204 207.710000, 120.218970 136.770000, and 147.208625 99.770000. The 'Search unmatched peaks only' checkbox is checked and circled in red. The 'Results' section has 'Overview' checked and 'Report top' set to '10 hits'. There are buttons for 'Copy Peaklist', 'Copy Masslist', 'Save as default', 'Start', and 'Exit'.

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

The image shows the Biotools software interface. On the left, the 'File' menu is open, and 'Combine multiple LIFT spectra...' is highlighted with a red circle. A green arrow points from this menu item to the 'Combine spectra' dialog box on the right. The dialog box has 'Input files' and 'Output file' sections. The 'OK' button is circled in red. A green arrow points from the 'OK' button to the 'Combine Spectra' text below.

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The image shows the Biotools software interface with a green arrow pointing to the 'PMF\_LIFT.mgf' file name in the bottom right corner of the window.

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